






AT-PMSC-I

LOW-TEMPERATURE WIDELINE AUTOTUNING
AND GONIOMETER NMR PROBE SYSTEM





Are you tired of manually tuning your NMR probe? Are you in possession of a precious sample containing exotic nuclei with spectral line widths of several MHz, broadened by electric quadrupolar or magnetic interactions? Does the recently discovered phase transformation in your frontier research multi-ferroic system take place at cryogenic temperatures? You want to study its exciting new physics using NMR, only to discover that your commercial probehead was not designed for low temperature applications? Well, if any of the above is true and you do not have engineering resources to design and manufacture a dedicated system yourself, then you may consider TEHO's new AT-PMSC-1 low-temperature wideline autotuning and goniometer NMR probe system.




See more at www.teho.com

ATC-400 Digital Autotuning Controller

Automated wideband NMR

AT-PMSC-1 is an automatic tuning and matching system for NMR spectroscopy, consisting of a wideband research-grade goniometer probe, detachable stepper motor drive unit, high/low power RF switch, and CompactPCI-based controller unit with touchscreen interface. It provides for automatic tuning of RF resonant circuits in NMR experiments with X-nuclei resonance frequencies up to 400 MHz. Its application areas span from wideband NMR of disordered solids where spectra are acquired using point-by-point frequency scan (frequency-stepping) methods, to investigations of smeared phase transitions via NMR of low-sensitivity nuclei or low-concentration samples where detuning takes place due to temperature-dependent capacitance and inductance of the trimming elements of the NMR detection circuit.



Remote control

If the autotune functionality is needed in the course of experiment, it can be controlled via the RS-232/Ethernet connection from the workstation which controls the operation of the spectrometer. Specifically, the synchronization of automatic tuning of AT-PMSC-1 with the spectrometer is established through AT-PMSC-1 procedure calls from within the spectrometer automation sequences (e.g. Python or C-based automation programming). This functionality is provided with all modern spectrometers.

Easy integration

AT-PMSC-1 is typically integrated into the RF detection circuit of the spectrometer by an in-line connection of the remotely controlled RF switch, immediately in front of the probe. The RF switch is used to redirect the RF line of the NMR probe from the autotuning controller to the NMR spectrometer during the signal acquisition mode. In this way, operation of AT-PMSC-1 that is independent of spectrometer's hardware is achieved; no RF or I/O connectivity between the spectrometer and AT-PMSC-1 is required.



Modular design

AT-PMSC-1 is an automatic tuning and matching system for NMR spectroscopy which consists of four components:

- ATP-CF wideband multipurpose probe, equipped with up to four trimmer capacitors (room-temperature range or cryogenic-temperature range versions) and high precision goniometer
- ATM-C4G rapidly detachable stepper motor unit for mounting onto ATP-CF, with trimmer status LED display panel
- ATS-500 RF switch with status LEDs
- ATC-400 digital controller unit with touchscreen display



The ATP-CF probe is designed to fit into the Oxford Instruments Spectrostat^{NMR} wide-bore cryostats of $\Phi = 2.5''$ (model ATP-63G-CF) or $\Phi = 2''$ (model ATP-51G-CF) bore inner-diameter. The probe is built to last: the housing material is highly corrosion-resistant non-magnetic stainless steel, the probehead base plate is made of titanium and the goniometer is machined from Vespel, delivering an exceptional wear resistance and hassle-free performance even at cryogenic temperatures. Stainless-steel ball-bearings are used to seat the trimmer and goniometer vacuum-sealed shafts. The mating shaft adap-

tors allow for a quick attachment or detachment of the ATM-C4G stepper motor unit. When the stepper motor unit is detached from the ATP-CF or not powered, the ATP-CF probe, specifically the trimmer and goniometer shafts, can be rotated manually for setup purposes. The ATM-C4G stepper motor unit is designed with industry-standard high power hybrid stepper motors and differential encoders so that it can be controlled with conventional stepper motion controllers. Your expensive, high-Q, cryogenic trimmer capacitors are protected from mechanical damage by adjustable limit switches.




Custom applications

Special care has been taken to ensure customizability of the probe, e.g. use of solenoids of various sizes and shapes, easy replacement of trimmer capacitors, and installation of additional mechanical, electrical, and optical components that may be needed in a custom-designed experiment. This is achieved via vacuum sealable access holes in the probe housing with pre-installed vacuum-tight multipole connectors with internal pinning, providing for easy attachment of RF circuits and wiring of external sensors, fiber-optic elements, and power supplies. The installation of user-designed coils has been made easy through the use of a rapid-exchange coil holders. The standard NMR detection circuit featuring an irradiation/detection coil and two trimming elements for tuning and matching can be upgraded to four trimming elements. Specifically, a second pair of variable capacitors can be installed, allowing for even broader frequency tuning range in the frequency-



stepping experiments (typically beyond 20 MHz) or for double resonance experiments. The use of solder has been reduced to a minimum when assembling the probe parts. Nevertheless, care has been taken to retain the required level of electromagnetic shielding performance. All mechanical parts can be disassembled using only a screwdriver. For a researcher with engineering skills, this makes it simple to alter the existing parts or to add custom-designed probe parts in order to expand the functionality of the probe.





per motor unit is powered directly from the ATC-400 via a single cable harness. The controller also provides power, control, and RF tuning signal for the ATC-500 RF switch via two additional cables. Complete wiring of the AT-PMSC-1 autotuning system into an existing NMR spectrometer is thus achieved by solely three, ready to use, flexible cables. Moreover, cabling lengths can be as much as 20m or more since optically-coupled and differential-type signalling is employed. Distant positioning of the NMR magnet, across a large lab from the NMR spectrometer, therefore presents no obstacle for the installation of the AT-PMSC-1 autotuning system.

Versatile operation

The ATC-400 controller unit can be used as a standalone device, i.e. as a wobbler, to tune any RF resonant circuit. This is particularly useful in the preparation phase of an NMR experiment whenever user designed RF coils are used in order to improve the S/N ratio or to expand the frequency range of the probe tunability. The controller unit is based on the industrial-grade CompactPCI platform fitted with Intel Core 2 Duo CPU blade and digital RF synthesizer. A separate internal power supply provides power to integrated RF reflectometer and stepper motor controller and drive units. The ATM-C4G step-






Touchscreen user interface

Through the ATC-400 touchscreen user interface, up to 4 trimmers and the goniometer of the ATP-CF probe can be controlled. The main screen of the interface displays a reflected RF power level vs. frequency and the orientation of the goniometer in real time, as well as the status bar. There are two modes of operation, the local one which, in addition to automatic tuning, allows for manual adjustment of trimmers and goniometer, and the remote one, intended for spectrometer-controlled autotuning and reorienting of samples. ATC-400 settings are controlled through a separate screen, typically used to select active trimmers, to initialize the positions of trimmer and goniometer shafts, and to access the

online help. Although the application software runs on top of Microsoft Windows XP, the OS is hidden from the operator's access. The complete ATC-400 software layer can be automatically restored or upgraded by inserting a bootable USB key when the system is powered-on. Frequency-stepping experiments in a broad frequency range (several MHz) as well as temperature-dependence experiments in a broad temperature range (several hundred K), which require precise tuning at each step, are easily performed by issuing string-based commands to the RS-232 bus from the automation macro programs being executed on the NMR spectrometer.



Technical specifications



ATC-400 digital controller

General

dimensions (W x D x H) without feet and handles	44.9 cm x 37.5 cm x 31 cm (17.7" x 14.8" x 12.2") industry-standard 19" 7U rack mountable, EMI shielded case
weight	25 kg (55 lbs)
housing material	chromated aluminium, surface paint in RAL 7038 color
power input rating	125/250 V AC, 50/60 Hz, 3 A

RF reflectometer

frequency sweep generation	DDS
frequency range	10 - 400 MHz
accuracy	1 ppm
output impedance	50 Ω
reflectometer range	0 - 50 dB

Stepper motor control

no. of channels	5 (4 x trimmer with limit switches, 1 x goniometer)
max. current	1.5 A/channel
inputs	2 x limit switch/channel, 1 x differential encoder/channel, motor status, emergency stop

Data acquisition and processing hardware

platform	32-bit CompactPCI
CPU unit	Intel Core Duo 2.26 GHz
DAQ	PCI bus, 16-bit, 1MS/s sample rate
remote connectivity	RS-232, Ethernet


Operator panel (front)

switches	on/off switch with two-color LED indicator
status indicators	standby mode LED remote mode LED
connectors	3 x USB 2.0
display panel	LCD, resistive touchscreen display size (diagonal): 12.1", 16:10 resolution: 1280 x 800

Operator panel (back)


switches	main power on/off switch
connectors	line power with fuse holder and line filter trimmers 1 - 4 : 19-pin female, color-coded goniometer: 19-pin male, color-coded ATS (to RF switch): 6-pin female ATM (to stepper motor unit): 4-pin female RF tune: type N female 1 x USB 2.0 1 x Ethernet 1 x RS-232

Software	
background OS (hidden from the operator)	Microsoft Windows XP SP3
main screen	hardware status line reflected power level vs. frequency plot panel with zoom buttons and level indicators input/display fields for tuning frequency (with lock) and goniometer orientation manual tuning and automatic tuning buttons mode selection buttons
settings screen	active trimmers selection panel trimmer drive mechanism initialization buttons goniometer orientation initialization button help system access button OS and remote configuration buttons
remote control	via RS-232 or Ethernet, string-based command language



ATM-C4G stepper motor unit

General	
max. diameter	Φ 13 cm (Φ 5.1")
full height (including supports)	34 cm (13.4")
weight	3.2 kg (7 lbs)
materials	housing: chromated aluminum, surface paint in RAL 7038 color supports: non-magnetic stainless steel limit switch compartment cover: transparent acrylic glass
power input rating (logic circuits)	+24 V, 500 mA DC
Stepper motors	
number of axes	3 - 5 (2 - 4 with low/high limit sensors, goniometer axis without limit sensors)
max. current	0.67 A per axis, powered from the autotuning controller
holding torque	127 Nmm (18 oz-in)
size and winding type	NEMA-11, bipolar serial
overheat protection	microcontroller-based with fan
Encoders	
resolution	500 steps/revolution
zero-position signal	yes
cable drivers	differential
Limit switch compartment	
limit switch type	mechanical microswitch with adjustable position
variable capacitor support	max. 100 turns full range
position indication	sliding limit switch trigger ± 50 turns, color-coded scales
Operator panel	
switches	on/off switch with LED indicator emergency switch
status indicators	motor activity LEDs (5 x) limit switch LEDs (8 x)
connectors	trimmers 1 - 4 : 19-pin female, color-coded goniometer: 19-pin male, color-coded power & status: 4-pin female



ATP-CF probe

General

compatibility	for use with Oxford Instruments Spectrostat ^{NMR} continuous flow CF-1200 cryostat
max. diameter	Φ 61.5 mm (Φ 2.42"), model ATP-63G-CF Φ 48.8 mm (Φ 1.92"), model ATP-51G-CF
full length (including shafts)	30 cm - 120 cm (11.8" to 47.2"), model ATP-63G-CF 30 cm - 60 cm (11.8" to 23.6"), model ATP-51G-CF exact dimensions are to be provided by the customer
weight	3.5 kg (7.7 lbs) at 90 cm (35.4") length, model ATP-63G-CF 2.5 kg (5.5 lbs) at 50 cm (19.7") length, model ATP-51G-CF
materials	housing: stainless steel, non-magnetic, high corrosion resistance goniometer: Vespel probehead base: titanium
temperature range	4K - 400 K, limited by the operating temperature range of trimmer capacitors
trimmer/goniometer adjustment	manual or stepper motor-controlled mechanical coupling to drive unit via male shaft couplers
vacuum sealing	rotary components (shafts): spring loaded bearing seals connectors: gasket seals

RF resonant circuit


type	parallel, 2 - 4 trimmer capacitors
trimmers	variable trimmer capacitors, standard or cryogenic version, extended shaft 0.5 pF min. capacitance, 70 pF max. capacitance, various capacitance ranges and working voltages available 20 - 105 turns full range max. shaft torque 70 Nmm (10 oz-in) 76 mm (3") max. length, Φ 18 mm (Φ 0.7") max. diameter
trimmer mounting	two grounded mounting holes for tuning capacitors two isolated (Teflon or Vespel) mounting holes for matching capacitors
trimmer rotational life	10,000 to 100,000 turns
irradiation/detection coil ^(*)	solenoid, saddle or other shapes rapid exchange coil holders (two pin male plug, Teflon base) for coil soldering Φ 20 mm (Φ 0.8") max. sample diameter 30 mm (1") max. sample length, model ATP-63G-CF 20 mm (0.8") max. sample length, model ATP-51G-CF
coil mounting	female socket, Teflon isolation, fixed into titanium probehead base
RF conductor	2 x semi-rigid, 50 Ω copper coaxial cables
external connectors	2 x N type female
internal connectors to trimmers	2 x SMA female, replaceable male counterpart with bare rigid cable termination

Goniometer

gear ratio	36:1 (10°/revolution)
resolution	0.02°
positioning accuracy	± 0.5°
sample holder seating	hexagonal, with fixing screw
max. sample holder dimensions	^(*) see coil data above

Auxiliary inputs

sensor connector S1	5-pin female
power connector S2	2-pin male, I _{max} = 4 A, V _{max} = 1 kV
spare hole	Φ 10 mm, vacuum-sealed when not used



ATS-500 high/low power RF switch

General

dimensions (W x D x H)	10.5 cm x 10.5 cm x 5 cm (4.1" x 4.1" x 2")
weight	0.5 kg (1.1 lbs)
housing material	chromated aluminium, surface paint in RAL 7038 color

Electrical characteristics

frequency range	DC – 2 GHz
max. CW power	400 W @ 500 MHz
max pulsed power	> 1.5 kW @ 500 MHz
impedance	50 Ω
actuator	electromechanical, failsafe, unused port terminated (5 W max.)
actuation signal	TTL
max. stray magnetic field	10 G (1 mT)





Mechanical characteristics

switching speed	< 20 ms
life	1.000.000 cycles

Operator panel

status indicators	high/low power LEDs
connectors	probe: type N female RF measure: type N female RF tune: type N female ATC (to autotuning controller): 6-pin female

Notes



KMZ Zalar Miran s.p.
Predjamska cesta 12
1000 Ljubljana, SLOVENIA

Jeklotehna TEHO d.o.o.
Cesta na Brdo 85
1000 Ljubljana, SLOVENIA

Tel. +386 1 423 54 50 office | +386 31 689 308 mobile
e-mail: nmr@teho.com | web: www.teho.com